

Automatic 24/7 Evaporation



Leading Safety Standards

Superior Ease of Use

Reduced Cost of Ownership


 **heidolph**
research made easy



Do you require a system to evaporate unlimited volumes in a continuous and unattended process 24 hours a day/7 days a week?

Automatic 24/7 Evaporation

Save time by distilling large volumes in a continuous and unattended manner overnight and during the weekend



General

- Automatic refill of the evaporation flask with the solvent mixture according to your settings. Continuous automatic discharge of condensate
- The automatic module is compatible with all bench-top and large-scale evaporators from Heidolph
- The worldwide unique automatic residue drainage offers unattended distillation for days
- The storage and the collecting tank can be replaced during the process

Leading Safety Standards

- The complete system is sensor-controlled and **turns off the automatic module as well as all connected equipment completely in case of an emergency**
- The system **can be placed under a fume hood** which ensures the highest safety level for the operator
- An easily accessible emergency-button allows you to **stop the process immediately** by the operator

Superior Ease of Use

- Take advantage of the wide range of possibilities to **adjust the settings according to your individual needs**
- The process can be **terminated automatically** after evaporation is completed or manually at any time by the operator
- Integrated plug connectors enable you to connect and to **control** also peripheral equipment **from a central point**
- The **worldwide unique automatic residue drainage** offers a wide range of new possibilities for **automatic, continuous and unattended distillation**

Reduced Cost of Ownership

- **Affordable alternative to common 50 or 100 liter evaporators**
- **Increase your efficiency** - the fully automatic process avoids the need of visual control and gives you time to focus on other tasks
- **Reduce your investments** by using your existing equipment - the automatic module is compatible to all rotary evaporators from Heidolph
- The unattended and continuous reconditioning of solvent **avoids any bottlenecks** of availability of solvents and reduce your investments in this way

➤ The smart way to process large quantities automatically

- The system is available for Hei-VAP bench-top evaporators and Hei-VAP Industrial evaporators
- Bumping is eliminated by continuously feeding small volumes
- For distillation select between these two modes:
 - **Sensor-Mode:** Refill and drainage are automatically sensor-controlled
 - **Time-Mode:** All parameters controlled by individually pre-programmed automatic time intervals
- Process rates with bench-top evaporators for the most common solvents are:
 - **Water: 1.4 l/h**
 - **Acetone: 6.0 l/h**
 - **Toluene: 6.5 l/h**
 - **Ethanol: 2.5 l/h**
- Increased process rates with Hei-VAP Industrial large-scale evaporators, for example:
 - **Water: 4.5 l/h**
 - **Acetone: 25.0 l/h**
- Both modes are completed in an automated and unattended process of solvent separations of nearly 100 %
- Degree of concentration for final product can be determined to rule out dehydration, crystallisation and formation of powder

➤ For distillation select between these two modes:



Sensor-Mode:

Refill and drainage are automatically sensor-controlled

- A sensor detects the fill-level inside the evaporation flask
- The desired/requested fill-level can be set individually
- Flask refill is achieved via valve control and under system pressure
- A sensor monitors the quantity of the solvent mixture and powers off the unit once it is completely processed

Time-Mode:

All parameters controlled by individually pre-programmed automatic time intervals

- The desired/requested fill-level can be set individually
- Flask refill is achieved via valve control and under system pressure according to pre-programmed cycles
- The process is powered off once the pre-programmed number of cycles have expired

Automated 24/7 Evaporation



Leading Safety Standards

Superior Ease of Use

Reduced Cost of Ownership

The average operational **lifespan of 10 years** is backed by a **3 year warranty** and makes your purchase a worthwhile investment.

Automated in unattended distillation in **two different modes**

Automatic refill of the evaporation flask and continuous **automatic discharge of condensate**

Automatic drainage of residue: residue in the evaporation flask will be drained off

Compatible with all Heidolph bench-top and large-scale rotary evaporators

The bench-top system can be **used in common fume hoods** due to the space saving design

Applications are completed in an automated and unattended process of **solvent separations of nearly 100 %**



Wide range of possibilities to **adjust settings according to your individual needs**

In case of an emergency **the complete system will power off entirely**

Bumping is eliminated by continuously feeding small volumes

➤ Functional Principle

For unattended and continuous use – including automatic drainage of residue

Save time by distilling large volumes in a continuous and unattended manner overnight and during the weekend

Automatic solvent supply of the evaporating flask

- Refilling with the Sensor-Mode runs automatically and is based on the desired level in the evaporation flask
- Refilling with the Time-Mode is done according to individual time parameters

Automatic discharge of condensate

- A valve matrix controls the continuous discharge of the condensate from the system into a temporary storage flask. Afterwards the collected condensate will be discharged into the collecting tank by a diaphragm pump

End of distillation process

- The process can be terminated automatically after evaporation is completed or manually at any time by the operator

Automatic drainage of residue

- The residue in the evaporating flask will be drained into a collecting tank. Afterwards the process starts again automatically by refilling the evaporating flask

➤ Distimatic Bench-top / Distimatic Industrial

The automatic module Hei-VAP Distimatic is compatible with all bench-top and large-scale evaporators from Heidolph



Distimatic Bench-top including automatic drainage of residue
P/N 591-53000-00

Distimatic Bench-top without automatic drainage of residue
P/N 591-51000-00



Distimatic Industrial including automatic drainage of residue and adaptor piece for glassware set A / A2
P/N 591-54000-00


Distimatic Industrial including automatic drainage of residue and adaptor piece for glassware set R
P/N 591-55000-00

Distimatic Industrial without automatic drainage of residue
P/N 591-52000-00

➤ Functional Principle

Distimatic Bench-top including automatic drainage of residue

- 1 Vacuum valve
- 2 Collector valve
- 3 Ventilation valve
- 4 Condensate pump
- 5 Residue pump



Rate of evaporation at continuous operation

- Reach outstanding evaporation rates of up to **6.5 l/h** with the **Distimatic Bench-top** and up to **25.0 l/h** with the **Distimatic Industrial**

› Distimatic Bench-top Packages

Distimatic Bench-top Silver Packages	P/N
Distimatic Bench-top Silver 1 Package Distimatic Bench-top including automatic drainage of residue, Hei-VAP Value Rotary Evaporator with coated, vertical glassware G9B	560-01914-00
Distimatic Bench-top Silver 2 Package Distimatic Bench-top including automatic drainage of residue, Hei-VAP Value Digital Rotary Evaporator with coated, vertical glassware set G9B	560-01914-20
Distimatic Bench-top Silver 3 Package Distimatic Bench-top without automatic drainage of residue, Hei-VAP Value Rotary Evaporator with coated, vertical glassware set G3B	560-01313-00
Distimatic Bench-top Silver 4 Package Distimatic Bench-top without automatic drainage of residue, Hei-VAP Value Digital Rotary Evaporator with coated, vertical glassware set G3B	560-01313-20
Distimatic Bench-top Gold Packages	P/N
Distimatic Bench-top Gold 1 Package Distimatic Bench-top including automatic drainage of residue, Hei-VAP Advantage Rotary Evaporator with hand lift and coated, vertical glassware set G9B	561-01914-00
Distimatic Bench-top Gold 2 Package Distimatic Bench-top including automatic drainage of residue, Hei-VAP Advantage Rotary Evaporator with motor lift and coated, vertical glassware set G9B	562-01914-00
Distimatic Bench-top Gold 3 Package Distimatic Bench-top without automatic drainage of residue, Hei-VAP Advantage Rotary Evaporator with hand lift and coated, vertical glassware set G3B	561-01313-00
Distimatic Bench-top Gold 4 Package Distimatic Bench-top without automatic drainage of residue, Hei-VAP Advantage Rotary Evaporator with motor lift and coated, vertical glassware set G3B	562-01313-00
Distimatic Bench-top Platinum Packages	P/N
Distimatic Bench-top Platinum 1 Package Distimatic Bench-top including automatic drainage of residue, Hei-VAP Precision Rotary Evaporator with hand lift and coated, vertical glassware set G9B	563-01914-00
Distimatic Bench-top Platinum 2 Package Distimatic Bench-top including automatic drainage of residue, Hei-VAP Precision Rotary Evaporator with motor lift and coated, vertical glassware set G9B	564-01914-00
Distimatic Bench-top Platinum 3 Package Distimatic Bench-top without automatic drainage of residue, Hei-VAP Precision Rotary Evaporator with hand lift and coated, vertical glassware set G3B	563-01313-00
Distimatic Bench-top Platinum 4 Package Distimatic Bench-top without automatic drainage of residue, Hei-VAP Precision Rotary Evaporator with motor lift and coated, vertical glassware set G3B	564-01313-00

› Distimatic Industrial Packages

Distimatic Industrial Platinum Packages	P/N
Distimatic Industrial Platinum 1 Package Distimatic Industrial including automatic drainage of residue and Hei-VAP Industrial with glassware set R	518-11155-00
Distimatic Industrial Platinum 2 Package Distimatic Industrial including automatic drainage of residue and Hei-VAP Industrial with glassware set A	518-13154-00
Distimatic Industrial Platinum 3 Package Distimatic Industrial including automatic drainage of residue and Hei-VAP Industrial with glassware set A2	518-14154-00
Distimatic Industrial Platinum 4 Package Distimatic Industrial without automatic drainage of residue and Hei-VAP Industrial with glassware set A2	518-14152-00
Distimatic Industrial Platinum 5 Package Distimatic Industrial without automatic drainage of residue and Hei-VAP Industrial with glassware set A	518-13152-00
Distimatic Industrial Platinum 6 Package Distimatic Industrial without automatic drainage of residue and Hei-VAP Industrial with glassware set R	518-11152-00

› Functional principle Distimatic



[Watch Video](#)

www.heidolph.de



➤ Technical Specifications - Automatic Module Distimatic

Model	Distimatic Bench-top without drainage of residue	Distimatic Bench-top including drainage of residue
P/N (230 V)	591-51000-00	591-53000-00
Compatibility (Rotary Evaporators)	Hei-VAP and Laborota 4000 series *	Hei-VAP and Laborota 4000 series
Separate control panel	yes	yes
Display	digital	digital
Automatic refill of sample	yes	yes
Automatic discharge of condensate	yes	yes
Automatic drainage of residue	no	yes
Power input control box (VA)	50	50
Noise level (dB(A))	< 80	< 80
Total power multiple switched socket (W)	3,500	3,500
Weight (kg)	5	5
Measurements without mounted components (W x D x H) (mm)	300 x 200 x 120	300 x 200 x 120
Permissible ambient temperature	5 – 31 °C at 80 % rel. humidity 32 – 40 °C linearly declining to max. 50 % rel. humidity 0 - 2000 m a.s.l., Contamination level 2 Overvoltage category II	5 – 31 °C at 80 % rel. humidity 32 – 40 °C linearly declining to max. 50 % rel. humidity 0 - 2000 m a.s.l., Contamination level 2 Overvoltage category II
Protection class (DIN EN 60529)	IP 54	IP 54

Model	Distimatic Industrial without drainage of residue	Distimatic Industrial including drainage of residue
P/N (230 V)	591-52000-00	591-54000-00 Glassware A / A2 591-55000-00 Glassware R
Compatibility (Rotary Evaporators)	Hei-VAP Industrial and Laborota 20 series	Hei-VAP Industrial and Laborota 20 series
Separate control panel	yes	yes
Display	digital	digital
Automatic refill of sample	yes	yes
Automatic discharge of condensate	yes	yes
Automatic drainage of residue	no	yes
Power input control box (VA)	50	50
Noise level (dB(A))	< 80	< 80
Total power multiple socket (W)	3,500	3,500
Weight (kg)	5	5
Measurements without mounted components (W x D x H) (mm)	300 x 200 x 120	300 x 200 x 120
Permissible ambient temperature	5 – 31 °C at 80 % rel. humidity 32 – 40 °C linearly declining to max. 50 % rel. humidity 0 - 2000 m a.s.l., Contamination level 2 Overvoltage category II	5 – 31 °C at 80 % rel. humidity 32 – 40 °C linearly declining to max. 50 % rel. humidity 0 - 2000 m a.s.l., Contamination level 2 Overvoltage category II
Protection class (DIN EN 60529)	IP 54	IP 54

All shown values relate on Automatic Module Disitmatic without any external accesory parts
* Also compatible to some models of other manufacturers. To check the compatibility please contact your dedicated Heidolph contact person


Certificate

To confirm the ability for
continuous operation
of the Automatic Module Distimatic

The Automatic Module Distimatic features overtemperature safety circuits according to DIN EN 61010-1:2001 and DIN EN 61010-2-010:2003 and therefore is designed for continuous operation.

This statement is made under the precondition that all units are operated in accordance with the operation manual and in accordance with good practice standards for safety in laboratories, rules for accident preventions, and compliance with directions on hazardous materials.

S c h w a b a c h , J a n u a r y 2 0 1 6


Stefan Peters
Research and Development Manager


Marcell Sarré
Quality Manager



Heidolph Instruments GmbH & Co. KG

Walpersdorfer Str. 12 • 91126 Schwabach/Germany

Phone +49 9122/9920-19 • Fax +49 9122/9920-65

Carina.Wagner@heidolph.de • www.heidolph.com

 **heidolph**
research made easy